



ALTAMEDICA INCORPORATED

Product List CHEMISTRY

Chemistry Analyzers

Code	Analyzers	Description	Brand
	Full Automated		
1805500	LIDA 500	420 T/hour, 540 T/hour with ISE; Sample positions: 100, Reagent positions: 80; 100 cuvettes with washing station; 2 reagent probes, 2 mixing probes; Light source: Halogen Tungsten Lamp	Linear Chemicals
1805000	LIDA 300	200 T/hour, 330 T/hour with ISE; Sample positions: 40, Reagent positions: 40; 40 cuvettes with washing station; 1 reagent/sample probe; Light Source: Halogen-Tungsten Lamp	Linear Chemicals
1800100	KROMA PLUS	240 T/hour, 360 T/hour with Ise; Sample positions: 49, Reagent positions: 35; 80 cuvettes with washing station; 1 reagent/sample probe (2 reagent/sample probe optional for upgrade); Light source: Halogen lamp	Linear Chemicals
1800050	KROMA	150 T/hour, 250 T/hour with ISE; Sample positions: 15, Reagent positions: 20; 80 cuvettes with washing station; 1 reagent/sample probe; Light source: Halogen lamp	Linear Chemicals
N/A	EXC200	240 T/hour for single reagent, 160 T/hour for dual reagent; Sample positions: 40, Reagent positions: 40; 63 cuvettes with washing station; 1 reagent/sample probe and 1 mixer probe; Light source: Halogen lamp	Zybio Inc.
N/A	CHEMRAY 240	240 T/hour; Sample positions: 40, Reagent positions: 40; 81 cuvettes, 9 cuvettes/stip with washing station; 1 reagent/sample probe, 1 mixing probe; Light source: Halogen Lamp	Rayto
N/A	CHEMRAY 120	100 T/hour; Sample positions: 9 (by default but can be adjusted up to 42 positions since the sample and reagent positions can be manually adjusted) Reagent positions: 26 (by default); 90 disposable cuvettes, 10 cuvettes/strip; 1 Reagent/sample probe, 1 mixing probe; Light source: Halogen lamp	Rayto
	Semi-Automated		
1803050	TEMIS	Wavelength range: 340nm-700nm w/ 2 free positions; 130 programs capacity; 600 patients results automatically stored; Light Source: Halogen lamp; With Internal thermal printer	Linear Chemicals



ALTAMEDICA INCORPORATED

Product List
CHEMISTRY

Code	Analyzers	Description	Brand
N/A	RT-9200	Wavelength range: 340nm-620nm, with 3 filters optional; 60 programs capacity; 2,200 patient results automatically stored; Light source: Halogen lamp; With Internal thermal printer	Rayto
N/A	RT-1904C	Wavelength range: 340nm-670nm with 1 filter optional; 160 programs capacity; 500 patient results automatically stored; Light source: Halogen lamp; With Internal thermal printer	Rayto
N/A	RT-9600	Wavelength range: 330nm-800nm with 1nm step-by adjustment; 200 programs capacity; 10,000 patient results automatically stored; Light source: Halogen lamp; With Internal thermal printer	Rayto
N/A	C61	Wavelength range: 340nm-620nm; 60 programs capacity; 2,000 patient results automatically stored; Light source: Halogen lamp; With Internal printer	Benesphera

Washing Solution

Code	Product Name	Packaging	Shelf-Life
1185010	LIDA Cleaner Solution	1 Liter	29 months
1185015	LIDA Extra-Wash Cuvette Solution	10x60 mL	29 months
KR18100	Kroma Extra-Wash Cuvettes Solution	8x32 mL	29 months
KR18300	Kroma Systemic Solution	10x50 mL	29 months
KR18400	Kroma Deproteinization Solution	1x500 mL	29 months
1180010	Iso-Clean	1x500 mL	29 months



ALTAMEDICA INCORPORATED

Product List
CHEMISTRY

Reagents

**Standard Included*

KR for Kroma Series and Lida 300

CT for Lida 500

Code	Reagent Name	Packaging	Availability
1100005	Acid Phosphatase	4x10 mL	Order Basis
1101000	Albumin*	2x50 mL	On Hand
1101010	Albumin*	4x100 mL	O.H.
1901005	Albumin Standard 5g/dL	1x3 mL	O.B.
KR10020	Albumin	2x30 mL	O.H.
KR10022	Albumin	8x30 mL	O.H.
CT10020	Albumin	2x40 mL	O.B.
1103005	Alkaline	2x50 mL	O.H.
1103010	Alkaline Phosphatase BR	3x100 mL	O.H.
KR10042	Alkaline Phosphatase BR	5x40 mL	O.B.
KR10045	Alkaline Phosphatase BR	11x40 mL	O.B.
CT10042	Alkaline Phosphatase BR	3x50 mL	O.B.
CT10045	Alkaline Phosphatase BR	12x50 mL	O.B.
1105000	ALT/GPT BR opt.	2x50 mL	O.H.
1105010	ALT/GPT BR opt.	3x100 mL	O.H.
KR10052	ALT/GPT BR opt.	5x40 mL	O.H.
KR10055	ALT/GPT BR opt.	11x40 mL	O.H.
CT10052	ALT/GPT BR opt.	3x50 mL	O.B.
CT10055	ALT/GPT BR opt.	12x50 mL	O.B.
1107005	Amylase MR	5x20 mL	O.H.
KR10060	Amylase MR	2x30 mL	O.B.
KR10062	Amylase MR	8x30 mL	O.B.
CT10060	Amylase MR	2x40 mL	O.B.
CT10062	Amylase MR	6x40 mL	O.B.
1109000	AST/GOT BR opt.	2x50 mL	O.H.
1109010	AST/GOT BR opt.	3x100 mL	O.H.
KR10072	AST/GOT BR opt.	5x40 mL	O.H.
KR10075	AST/GOT BR opt.	11x40 mL	O.H.



ALTAMEDICA INCORPORATED

Product List

CHEMISTRY

Code	Reagent Name	Packaging	Availability
CT10072	AST/GOT BR opt.	3x50 mL	O.B.
CT10075	AST/GOT BR opt.	12x50 mL	O.B.
1110005	Bilirubin Direct	2x50 mL	O.B.
1111010	Bilirubin Total	2x100 mL	O.B.
1112005	Bilirubin Total and Direct	2x100 mL	O.H.
1912005	Bilirubin Standard (T) (D)	1x1 mL	O.H.
KR10082	Bilirubin Direct	5x40 mL	O.B.
KR10092	Bilirubin Total	5x40 mL	O.B.
KR10095	Bilirubin Total	11x40 mL	O.B.
CT10082	Bilirubin Direct	3x50 mL	O.B.
CT10085	Bilirubin Direct	12x50 mL	O.B.
CT10092	Bilirubin Total	3x50 mL	O.B.
CT10095	Bilirubin Total	12x50 mL	O.B.
1113000	Calcium Arsenazo III Color*	2x50 mL	O.B.
1944005	Magnesium/ Calcium Standard	1x3 mL	O.B.
KR10110	Calcium Arsenazo III	2x30 mL	O.B.
KR10112	Calcium Arsenazo III	8x30 mL	O.B.
CT10110	Calcium Arsenazo III	2x40 mL	O.B.
CT10112	Calcium Arsenazo III	6x40 mL	O.B.
1115000	Calcium OCC*	2x50 mL	O.B.
1115010	Calcium OCC*	4x100 mL	O.B.
1944005	Magnesium/ Calcium Standard	1x3 mL	O.B.
1116005	Chloride*	2x50 mL	O.B.
1916005	Chloride/ Phosphorus Standard	1x3 mL	O.B.
KR10130	Chloride	2x30 mL	O.B.
CT10130	Chloride MR	2x40 mL	O.B.



ALTAMEDICA INCORPORATED

Product List
CHEMISTRY

Code	Reagent Name	Packaging	Availability
1118005	Cholesterol MR*	2x50 mL	O.H.
1118010	Cholesterol MR*	4x100 mL	O.H.
1118015	Cholesterol MR*	4x250 mL	O.H.
1918005	Cholesterol Standard 200mg/dL	1x3 mL	O.B.
KR10140	Cholesterol MR	2x30 mL	O.B.
KR10142	Cholesterol MR	8x30 mL	O.H.
KR10145	Cholesterol MR	18x30 mL	O.H.
CT10140	Cholesterol MR	2x40 mL	O.B.
CT10142	Cholesterol MR	6x40 mL	O.B.
CT10145	Cholesterol MR	10x60 mL	O.B.
1119005	Cholinesterase (Total and Inhibited)	30 Tests	O.B.
1120010	Creatinine Kinase BR	1x25 mL	O.B.
1120005	Creatinine Kinase BR	2x50 mL	O.B.
1923005	Creatinine Standard 2mg/dL	1x3 mL	O.B.
KR10160	Creatinine Kinase BR	1x40 mL	O.B.
KR10162	Creatinine Kinase BR	5x40 mL	O.B.
CT10160	Creatinine Kinase BR	1x50 mL	O.B.
CT10162	Creatinine Kinase BR	3x50 mL	O.B.
1123005	Creatinine Kinase MB	1x25 mL	O.B.
KR10170	Creatinine Kinase MB	1x40 mL	O.B.
KR10172	Creatinine Kinase MB	5x40 mL	O.B.
CT10170	Creatinine Kinase MB	1x50 mL	O.B.
CT10172	Creatinine Kinase MB	3x50 mL	O.B.



ALTAMEDICA INCORPORATED

Product List
CHEMISTRY

Code	Reagent Name	Packaging	Availability
1123005	Creatinine*	2x50 mL	O.H.
1123010	Creatinine*	4x100 mL	O.H.
1123020	Creatinine*	4x250 mL	O.H.
1923005	Standard 2mg/dL	1x3 mL	O.B.
KR10150	Creatinine E	5x45 mL	O.H.
KR10152	Creatinine E	11x45 mL	O.H.
CT10152	Creatinine E	3x45 mL	O.B.
CT10155	Creatinine E	12x45 mL	O.B.
KR10460	Fructosamine MR	2x30 mL	O.B.
CT19460	Fructosamine Standard	2x3 mL	O.B.
1126005	GGT BR opt.	2x50 mL	O.B.
1126010	GGT BR opt.	3x100 mL	O.B.
KR10190	GGT BR opt.	1x40 mL	O.B.
KR10192	GGT BR opt.	5x40 mL	O.B.
CT10192	GGT BR opt.	3x50 mL	O.B.
CT10195	GGT BR opt.	12x50 mL	O.B.
1129005	Glucose MR*	2x50 mL	O.H.
1129010	Glucose MR*	4x100 mL	O.H.
1129015	Glucose MR*	4x250 mL	O.H.
1929005	Glucose Standard 100mg/dL	1x3 mL	O.B.
KR10200	Glucose MR	2x30 mL	O.B.
KR10202	Glucose MR	8x30 mL	O.H.
KR10205	Glucose MR	18x30 mL	O.H.
CT10200	Glucose MR	2x40 mL	O.B.
CT10202	Glucose MR	6x40 mL	O.B.
CT10205	Glucose MR	10x60 mL	O.B.
1130010	GOT Color	2x100 mL	O.B.
1132010	GPT Color	2x100 mL	O.B.
1133505	HDL-Cholesterol Direct	40 mL	O.B.
1133510	HDL-Cholesterol Direct	320 mL	O.B.



ALTAMEDICA INCORPORATED

Product List
CHEMISTRY

Code	Reagent Name	Packaging	Availability
1133525	HDL-Cholesterol Direct	80 mL	O.H.
1972995	Human LDL/HDLc Calibrator	1x1 mL	O.B.
1972010	Human LPD Control	3x3mL	O.B.
KR10240	HDL-Cholesterol Direct	1x40 mL	O.B.
KR10242	HDL-Cholesterol Direct	5x40 mL	O.H.
KR10245	HDL-Cholesterol Direct	11x40 mL	O.H.
CT10240	HDL-Cholesterol Direct	1x60 ml	O.B.
CT10242	HDL-Cholesterol Direct	3x60 mL	O.B.
1133010	HDL-Cholesterol*	2x40 mL	O.H.
1933005	HDL-Cholesterol Standard 50mg/dL	1x3 mL	O.B.
1134015	Hemoglobin*	2x40 mL	O.B.
19334005	Hemoglobin Standard 12g/dL	1x1 mL	O.B.
	Hemoglobin Special	4x5 mL	O.B.
1163005	Hemoglobin MR	4x100 mL	O.B.
1934005	Hemoglobin Standard 12g/dL	1x1 mL	O.B.
1135105	Iron Cromazoul*	2x50 mL	O.B.
1938005	Iron Standard 100 ug/dL	1x3 mL	O.B.
CT10272	Iron Cromazoul	3x60 mL	O.B.
1135005	Iron Ferrozine*	2x50 mL	O.B.
1938005	Iron Standard 100ug/dL	1x3 mL	O.B.
KR10262	Iron Ferrozine	5x40 mL	O.B.
KR10265	Iron Ferrozine	11x40 mL	O.B.
CT10262	Iron Ferrozine	3x50 mL	O.B.
CT10265	Iron Ferrozine	12x50 mL	O.B.



ALTAMEDICA INCORPORATED

Product List
CHEMISTRY

Code	Reagent Name	Packaging	Availability
1141010	LDH BR	2x50 mL	O.B.
KR10282	LDH BR	5x40 mL	O.B.
CT10282	LDH BR	3x50 mL	O.B.
1142005	LDL-Cholesterol Direct	40 mL	O.B.
1142010	LDL-Cholesterol Direct	320 mL	O.B.
1972005	Human LDL/HDLc Calibrator	1x1 mL	O.B.
1972010	Human LPD Control	3x3 mL	O.B.
KR10290	LDL-Cholesterol Direct	5x40 mL	O.B.
KR10292	LDL-Cholesterol Direct	1x40 mL	O.B.
CT10290	LDL-Cholesterol Direct	1x60 mL	O.B.
1133105	LDL-Cholesterol*	1x5 mL	O.B.
1933105	LDL-Cholesterol	1x3 mL	O.B.
1143010	Lipase	1x60 mL	O.B.
KR10430	Lipase BR	1x36 mL	O.B.
1943005	Lipase Calibrator	1x3 mL	O.B.
1144005	Magnesium*	2x50 mL	O.B.
1944005	Magnesium/ Calcium Standard	1x3 mL	O.B.
KR10310	Magnesium BR	1x40 mL	O.B.
CT10310	Magnesium BR	2x40 mL	O.B.
CT10312	Magnesium BR	3x80 mL	O.B.
1149005	Phosphorus UV*	2x50 mL	O.B.
1149010	Phosphorus UV*	4x100 mL	O.B.
19166005	Chloride/ Phosphorus Standard	1x3 mL	O.B.
KR10320	Phosphorus UV	8x30 mL	O.B.
KR10322	Phosphorus UV	2x40 mL	O.B.
CT10320	Phosphorus UV	2x40 mL	O.B.
CT10325	Phosphorus UV	10x60 mL	O.B.



ALTAMEDICA INCORPORATED

Product List
CHEMISTRY

Code	Reagent Name	Packaging	Availability
CT10322	Phosphorus UV	6x40 mL	O.B.
1148010	Phosphorus*	3x50 mL	O.B.
1916005	Phosphorus Standard	1x3 mL	O.B.
1150015	Potassium*	1x25 mL	O.B.
1150010	Potassium*	3x50 mL	O.B.
KR10330	Potassium BR	1x40 mL	O.B.
CT10330	Potassium BR	1x50 mL	O.B.
CT10332	Potassium BR	3x50 mL	O.B.
1162005	Protein (Urine and CSF)	2x50 mL	O.B.
1162025	Protein (Urine and CSF)	1x250 mL	O.B.
1962005	Protein (Urine and CSF) Standard 2g/L	1x3 mL	O.B.
1151015	Sodium*	1x30 mL	O.B.
1151010	Sodium*	2x60 mL	O.B.
KR10340	Sodium BR	1x40 mL	O.B.
CT10340	Sodium BR	1x30 mL	O.B.
CT10342	Sodium BR	3x30 mL	O.B.
1137010	TIBC	50 Tests	O.B.
1153005	Total Protein*	2x50 mL	O.H.
1153010	Total Protein*	4x100 mL	O.H.
1153020	Total Protein*	4x250 mL	O.B.
1953005	Total Protein Standard 7g/dL	1x3 mL	O.B.
KR10350	Total Protein	2x30 mL	O.H.
KR10352	Total Protein	8x30 mL	O.H.
KR10355	Total Protein	18x30 mL	O.H.
CT10350	Total Protein	2x40 mL	O.B.
CT10352	Total Protein	6x40 mL	O.B.
CT10355	Total Protein	10x60 mL	O.B.
1155005	Triglycerides MR*	2x50 mL	O.H.



ALTAMEDICA INCORPORATED

Product List
CHEMISTRY

Code	Reagent Name	Packaging	Availability
1155010	Triglycerides MR*	4x100 mL	O.H.
1955005	Triglycerides Standard 200 mg/dL	1x3 mL	O.B.
KR10360	Triglycerides MR	2x30 mL	O.B.
KR10362	Triglycerides MR	8x30 mL	O.H.
KR10365	Triglycerides MR	18x30 mL	O.H.
CT10360	Triglycerides MR	2x40 mL	O.B.
CT10362	Triglycerides MR	6x40 mL	O.B.
CT10365	Triglycerides MR	10x60 mL	O.B.
1156010	Urea Berthelot*	2x50 mL	O.B.
1156015	Urea Berthelot*	4x100 mL	O.B.
1958005	Urea Standard 50 mg/dL	1x3 mL	O.B.
1158005	Urea/BUN BR	2x50 mL	O.H.
1158010	Urea/BUN BR	3x100 mL	O.H.
1958005	Urea Standard 50 mg/dL	1x3 mL	O.B.
KR10382	Urea/BUN BR	5x40 mL	O.H.
KR10385	Urea/BUN BR	11x40 mL	O.H.
CT10382	Urea/BUN BR	3x50 mL	O.B.
CT10385	Urea/BUN BR	12x50 mL	O.B.
1161005	Uric Acid MR*	2x50 mL	O.H.
1161010	Uric Acid MR*	4x100 mL	O.H.
1161015	Uric Acid MR*	4x250 mL	O.H.
1961005	Uric Acid Standard 6mg/dL	1x3 mL	O.B.
KR10390	Uric Acid MR	2x30 mL	O.B.
KR10392	Uric Acid MR	8x30 mL	O.H.
KR10395	Uric Acid MR	18x30 mL	O.H.
CT10390	Uric Acid MR	2x40 mL	O.B.
CT10392	Uric Acid MR	6x40 mL	O.B.
CT10395	Uric Acid MR	10x60 mL	O.B.

Disclaimer: All the products listed herein are subject to change from time to time without prior notice. Some items are pre-order basis. For inquiries concerning availability, packaging, and prices, kindly contact customerservice@altamedica.com.ph



ALTAMEDICA INCORPORATED

Product List CHEMISTRY

Chromatography in Tube

Code	Reagent Name	Packaging	Availability
3155105	Glycated HbA1c	25 Tests	Order Basis
3155110	Glycated HbA1c	26 Tests	Order Basis

Controls and Calibrators

Code	Control/ Calibrator	Packaging	Availability
1975005	Human Multicalibrator CC/H	5x5 mL	On Hand
1980005	Human Multisera Control Normal	5x5 mL	On Hand
1986005	Human Multisera Control Abnormal	5x5 mL	On Hand
CT19750	Human Multicalibrator CC/H	8x3 mL	Order Basis
CT19800	Human Multisera Control Normal	8x5 mL	Order Basis
CT19850	Human Multisera Control Abnormal	8x5 mL	Order Basis

Chemray 120 Consumables and Accessories

Code	Product Name	Packaging	Brand
1-3-070-048-XX	40ml Reagent Bottle	50 pcs/pack	Rayto
1-3-070-050-XX	40ml Reagent Bottle Cap (White)	50 pcs/pack	Rayto
1-3-070-047-XX	20ml Reagent Bottle	50 pcs/pack	Rayto
1-3-070-049-XX	20ml Reagent Bottle Cap (Blue)	50 pcs/pack	Rayto
1-3-080-005-XX	Test cuvette (1strip/10 rxn)	100strips/box	Rayto
N/A	Small Tube/cup	500pcs/pack	Rayto
N/A	RBA Cleanser	1L / bot	Rayto

Disclaimer: All the products listed herein are subject to change from time to time without prior notice. Some items are pre-order basis. For inquiries concerning availability, packaging, and prices, kindly contact customerservice@altamedica.com.ph



ALTAMEDICA INCORPORATED

Product List
CHEMISTRY

EXC200 Consumables and Accessories

Code	Product Name	Packaging	Brand
02-24-00-0145-00	Quartz Cuvette (6.6*6.6*32mm)	64pcs/pack	Zybio
02-65-00-0097-00	Quartz Cuvette (6.6*6.6*32 mm)	1pc	Zybio
02-24-00-0142-00	Plastic Cuvette	300pcs/pack	Zybio
02-36-10-0009-00	Bar code scanner	1pc	Zybio
02-24-00-0001-00	Light Source Assembly	1pc	Zybio
02-25-00-0209-00	Reagent Bottle, 20mL	100pcs/pack	Zybio
02-25-00-0208-00	Reagent Bottle, 35mL	100pcs/pack	Zybio
02-24-00-0146-00	Concentrated Detergent	5Lx4/box	Zybio
01.09.50.01.19.02	Probe Detergent	30mLx8/box	Zybio
02-24-00-0132-00	Waste Container & Pure Water Container	2pcs/pack	Zybio

Disclaimer: All the products listed herein are subject to change from time to time without prior notice. Some items are pre-order basis. For inquiries concerning availability, packaging , and prices, kindly contact customerservice@altamedica.com.ph