



## ALTAMEDICA INCORPORATED

### Product List HEMATOLOGY

#### *Hematology Analyzer*

Code	Product Name	Description	Brand
3802000	STEL 3	3-part differential hematology analyzer with 60 tests per hour throughput	Linear Chemicals
3803000	STEL 5	5-part differential hematology analyzer with 60 tests per hour throughput	Linear Chemicals
3803200	STEL5 Plus Autoloader Hematology Analyzer	5-part differential hematology analyzer with 80 tests per hour throughput	Linear Chemicals
N/A	Z3	3-part differential hematology analyzer with 70 tests per hour throughput	Zybio Inc.
N/A	Z3 with CRP	3-part differential hematology analyzer with 70 tests per hour for WB and 60 tests per hour for CRP throughput	Zybio Inc.
N/A	Hemaray 83	5-part differential hematology analyzer with 60 tests per hour throughput	Rayto
N/A	Hemaray 86	5-part differential hematology analyzer with 60 tests per hour throughput ( with autoloader)	Rayto

#### *Reagents and Controls for STEL 3*

Code	Product Name	Packaging	Shelf-Life
3802005	STEL3 Diluent	1 x 20L	2 years
3802010	STEL 3 Lyse	1 x 0.5L	2 years
3802015	STEL 3 Cleaner	1 x 1L	2 years
3802020	STEL 3 Extra Cleaner	1 x 0.1L	2 years
3802025	STEL 3 Control	3 (Levels) x 2.5mL	6 months

#### *Reagents and Control for STEL 5*

Code	Product Name	Packaging	Shelf-Life
3803005	STEL 5 Diluent	1 x 20L	2 years
3803010	STEL 5 LH Lyser	1 x 0.2L	2 years
3803012	STEL 5 Diff Lyse	1 x 0.5L	2 years
3803020	STEL 5 Extra Cleaner	1 x 50mL	2 years
3803025	STEL 5 Hematology Control	3(Levels) x 3mL	6 months



## ALTAMEDICA INCORPORATED

### Product List HEMATOLOGY

#### *Reagents, Control and Calibrator for STEL 5 Plus*

Code	Product Name	Packaging	Shelf-Life
3803210	STEL 5 Plus LD Lyser	1 x 1L	18 - 20 months
3803212	STEL 5 Plus LH Lyser	1 x 0.2L	18 - 20 months
3803205	STEL 5 Plus Diluent	1x 20L	18 - 20 months
3803220	STEL 5 Plus Probe Cleanser	1x 50mL	18 - 20 months
3803225	STEL 5 Plus Controls	3 x 3mL	6 months
3803230	STEL 5 Plus Cal	1 x 3mL	6 months

#### *Reagents and Control for Z3, Z3 with CRP, and Z5*

Code	Product Name	Packaging	Shelf-Life
N/A	Z3 Diluent	1 x 20L	2 years
N/A	Z3 Lyse	1 x 0.5L	2 years
N/A	Z3 Probe Cleanser	2 x 0.05L	1 year
N/A	CBC 3D for Zybion (Control for Z3)	3 (Levels) x 3mL	105 days
N/A	CBC CAL-PLUS for Zybion (Calibrator for Z3 and Z5)	1 x 3mL	45 days
N/A	CRP Reagent 1	25 tests per bag	1 year
N/A	CRP Reagent 2	2 x 3mL	1 year
N/A	CRP Control	2 (Levels) x 1mL	1 year
N/A	CRP Calibrator	6 (Levels) x 1mL	1 year
N/A	Z5 Diluent	1 x 20L	2 years
N/A	Z5 LD Lyse	4 x 0.5L	2 years
N/A	Z5 LB Lyse	4 x 0.5L	2 years
N/A	Z5 Probe Cleanser	2 X 50mL	1 year
N/A	CBC-DH for Zybion (Control for Z5)	3 (Levels) x 3mL	75 days
N/A	Z3 + Z5 Control Bundle	6 (3 levels) x 3mL	

#### *Hemaray83 Reagents*

Code	Product Name	Packaging	Brand
RD	Diluent	20L/box	Rayto
RL-H	Lyse H	0.5L/box	Rayto
RL-D	Lyse D	0.5L/box	Rayto
RC	Cleanser	1L/box	Rayto
REC	Concentrated Cleanser	0.03L	Rayto

#### *Hemaray86 Reagents*

Code	Product Name	Packaging	Brand
RD-86	Diluent	20L/box	Rayto
RT-86H	5H Lyse	0.5L/box	Rayto
RL-86D	5D Lyse	0.5L/box	Rayto
RC-86	Cleanser	1L/box	Rayto
RC-86C	Concentrated Cleanser	0.05L	Rayto



## ALTAMEDICA INCORPORATED

### Product List HEMATOLOGY

#### Controls and Calibrators

##### 3-Part Hematology Analyzer

Code	Product Name	Packaging	Brand
3D507	CBC-3D Tubes Hema Control (4L, 4N, 4H)	12x2mL	R & D
3D508	CBC-3D Tubes Hema Control (2L, 2N, 2H)	6x2mL	R & D
3D508	CBC-3D Tubes Hema Control (2L, 2N, 2H)	3x2mL	R & D

##### 5-Part Hematology Analyzer

Code	Product Name	Packaging	Brand
5D003	CBC-5D Tubes Hema Control 4L, 4N, 4H)	12x5mL	R & D
5D004	CBC-5D Tubes Hema Control 2L, 2N, 2H)	6x5mL	R & D
5D004	CBC-5D Tubes Hema Control 2L, 2N, 2H)	3x5mL	R & D
8CP33	Calibration	3x4.5mL	R & D

Disclaimer: All hematology controls are pre-order basis. For prices and schedule of ordering, kindly contact customerservice@altamedica.com.ph

#### Erythrocyte Sedimentation Rate (ESR) Machine

Code	Product Name	Description	Brand
5810005	Iria	ESR Machine with 5 position and a throughput of 18 tests per hour	Linear Chemicals
5810010	Lena 16 <sup>EDTA</sup>	ESR Machine with a throughput of 48 tests per hour and utilizes EDTA tube	Linear Chemicals
5810020	Lena 32 <sup>EDTA</sup>	ESR Machine with a throughput of 96 tests per hour and utilizes EDTA tube	Linear Chemicals
5810030	Lena NE	ESR Machine with a throughput of 50 tests per hour and utilizes Plus-SED or VacuSED tubes	Linear Chemicals
5810025	Therma NE	ESR Machine with a throughput of 250 tests per hour and utilizes Plus-SED or VacuSED tubes	Linear Chemicals
5310010	Visual-SED ESR Rack	ESR MANUAL RACK takes the analysis in 30 min for 1 hour results. Results are expressed in Westergren millimeters and it is capable to process 20 samples per hour.	Linear Chemicals



## ALTAMEDICA INCORPORATED

### Product List HEMATOLOGY

#### *Accessories for Erythrocyte Sedimentation Rate (ESR) Machine*

Code	Product Name	Packaging	Brand
N/A	Barcode Reader	1 piece	Linear Chemicals
N/A	Thermal Paper (57 x 250 cm)	1 roll	Linear Chemicals

#### *Controls and Accessories for Erythrocyte Sedimentation Rate (ESR) Machine*

Code	Product Name	Packaging	Shelf-Life
5910005	ESR Control (N-I N-II)	2 (Levels) x 10mL	
5920005	QC1: ESR Control (Plus-SED)	2 (Levels) x 1.5mL	
59310010	QC2: ESR (VacuSED)	2 (Levels) x 1.5mL	
5310020	Bar Code Reader		
5500240	Thermal Printer Paper	57x250cm	
5810092	LENA EDTA Dongle 2.000 Tests		
5220060	Plus-SED auto (Plastic tubes without Vacuum 1.5mL)	6 x 100 pcs.	
5231050	VacuSED (Plastic tubes with Vacuum 1.5mL)	5 x 100 pcs.	

#### *Coagulation Machines*

Code	Product Name	Description	Brand
3801000	Ares	A semi-automated coagulation analyzer with four channels.	Linear Chemicals
3801100	Ceres	A semi-automated coagulation analyzer with two channels.	Linear Chemicals
3801050	Althea	A fully-automated coagulation analyzer with a throughput of 60 tests per hour for PT and 50 tests per hour for PT and APTT	Linear Chemicals

#### *Consumables for Ares, Ceres, and Althea*

Code	Product Name	Packaging	Shelf-Life
3801920	Sample cups for Ares and Ceres	1 x 50 pcs.	N/A
3801820	Reaction cups for Althea	1 x 1000 pcs.	N/A
3801505	Althea Cleanser Solution	4 x 30 mL	3 years



## ALTAMEDICA INCORPORATED

### Product List HEMATOLOGY

#### *Reagents and Controls for Coagulation*

<b>Code</b>	<b>Product Name</b>	<b>Packaging</b>	<b>Shelf-Life</b>
3145005	D-Dimer Turbidimetric	1 x 54mL	18 months
KR31450	D-Dimer Turbidimetric	1 x 36mL	18 months
3945005	Calibrator for D-Dimer	2 x 1mL	18 months
3945010	Control for D-Dimer Turbidimetric	2 (Levels) x 0.5mL	18 months
3510101	PT-HS	10 x 5mL	2 years
3510205	PT-HS	5 x 5mL	2 years
3510201	APTT (Ellagic Acid)	10 x 4mL	2 years
3510205	APTT (Ellagic Acid)	5 x 4mL	2 years
3510301	Fibrinogen	5 x 2mL	2 years
3510401	Calcium Chloride 0.020M	3 x 16mL	2 years
3520101	Plasma Control Level 1	10x1mL	18 months
3520105	Plasma Control Level 1	5x1mL	18 months
3520201	Plasma Control Level 2	10x1mL	18 months
3520205	Plasma Control Level 2	5x1mL	18 months

*Disclaimer: All the products listed herein are subject to change from time to time without prior notice. For prices and availability, kindly contact customerservice@altamedica.com.ph*