



ALTAMEDICA INCORPORATED

Product List  
VETERINARY

*Rapid Tests*

*Canine Range*

Code	Product Name	Sample Type	Packaging	Shelf-Life
VE42011	CHW Ag (Heartworm)	WB/Serum/Plasma	10 tests per box	
VE42021	C.Bru-LPS (Brucella)	WB/Serum/Plasma	10 tests per box	
VE42031	LSH Ab (Leishmania)	WB/Serum/Plasma	10 tests per box	
VE42041	Ehrlichia	WB/Serum/Plasma	10 tests per box	
VE42051	cCRP (CRP)	WB/Serum/Plasma	10 tests per box	
VE42061	Ehrlichia Ab+LSH Ab	WB/Serum/Plasma	10 tests per box	
VE42071	Ehrlichia Ab+CHW Ag	WB/Serum/Plasma	10 tests per box	
VE42081	Ehrlichia Ab+LSH Ab+CHW Ag	WB/Serum/Plasma	10 tests per box	
VE42091	CPV Ag (Parvovirus)	Feces/Vomit	10 tests per box	
VE42101	CCV Ag (Coronavirus)	Feces/Vomit	10 tests per box	
VE42111	CPV+CCV Ag	Feces/Vomit	10 tests per box	
VE42121	CPV Ag+CCV Ag+Giardia Ag	Feces/Vomit	10 tests per box	
VE42131	CDV Ag (Distemper)	Secretions	10 tests per box	
VE42141	CAV Ag (Adenovirus II)	Secretions	10 tests per box	
VE42151	CIV Ag (Influenza)	Secretions	10 tests per box	
VE42161	CDV+CPV Ag	Secretions+ Feces/Vomit	10 tests per box	
VE42171	CDV Ab (Distemper)	Serum/Plasma	10 tests per box	

*Feline Range*

Code	Product Name	Sample Type	Packaging	Shelf-Life
VE43011	FeLV Ag (Leukemia)	WB/Serum/Plasma	10 tests per box	
VE43021	FIV Ab (Immunodeficiency)	WB/Serum/Plasma	10 tests per box	
VE43031	FIV Ab+FeLV Ag	WB/Serum/Plasma	10 tests per box	
VE43041	fSAA Ag (Amyloid A)	WB/Serum/Plasma	10 tests per box	
VE43051	FPV Ag (Panleucopenia)	Feces/Vomit	10 tests per box	
VE43061	FCV Ag (Calicivirus)	Secretions	10 tests per box	
VE43071	FIPV (Peritonitis)	Serum/Feces/ Ascites/Secretions	10 tests per box	
VE43081	FCoV (Coronavirus)	Pleural Fluid/Ascites	10 tests per box	
VE43091	FHV Ag (Herpes)	Eye Secretion	10 tests per box	



ALTAMEDICA INCORPORATED

**Product List**  
**VETERINARY**

*Porcine Range*

<b>Code</b>	<b>Product Name</b>	<b>Sample Type</b>	<b>Packaging</b>	<b>Shelf-Life</b>
VE44011	TRI (Trichinellosis)	WB/Serum/Plasma	10 tests per box	
VE44021	CSFV (Swine Fever)	Serum/Plasma	10 tests per box	
VE44031	PRRSV (Virus Ab)	Serum/Plasma	10 tests per box	
VE44041	P. FMDV-A (Virus A Ab)	Serum/Plasma	10 tests per box	
VE44051	PCV-2 (Circovirus)	Serum/Plasma	10 tests per box	
VE44061	PRV-gB (Pseudorabies) gB Ab	Serum/Plasma	10 tests per box	
VE44071	PRV-gE (Pseudorabies gE Ab)	Serum/Plasma	10 tests per box	

*Others Range*

<b>Code</b>	<b>Product Name</b>	<b>Sample Type</b>	<b>Packaging</b>	<b>Shelf-Life</b>
VE45011	Brucella Ag (Bovine)	WB/Serum/Plasma	10 tests per box	
VE45021	Toxo IgG IgM	WB/Serum/Plasma	10 tests per box	
VE45031	Giardia Ag	Feces	10 tests per box	
VE45041	Rota Antigen	Feces	10 tests per box	
VE45051	Toxo Ag	Feces/Serum/Plasma	10 tests per box	
VE45061	NDV Ag (Newcastle Disease Virus)	Feces/Secretions	10 tests per box	
VE45071	MPB 70&83 (Bovine Tuberculosis)	Serum/Plasma	10 tests per box	
VE45081	B. FMDV-A (Bovine Foot and Mouth Disease Virus A Ab Test)	Serum/Plasma	10 tests per box	
VE45091	BVDV (Bovine Viral Diarrhea)	Serum/Plasma	10 tests per box	
VE45101	CE Ab Test	Serum/Plasma	10 tests per box	
VE45111	FMDV-O (Foot and Mouth Disease Virus O Ab)	Serum/Plasma	10 tests per box	
VE45121	FMD NSP Ab	Serum/Plasma	10 tests per box	

*Urinalysis*

<b>Code</b>	<b>Product Name</b>	<b>Description</b>	<b>Packaging</b>	<b>Shelf-Life</b>
VE71001	Urinalysis Reagent Strips	LEU/NIT/URO/PRO/ pH/BLO/SG/KET/BIL / GLU/ASC	100 tests	
VE78000	Urine Analyzer			



ALTAMEDICA INCORPORATED

**Product List  
VETERINARY**

*Liquid Clinical Chemistry*

<b>Code</b>	<b>Product Name</b>	<b>Sample Type</b>	<b>Packaging</b>	<b>Shelf-Life</b>
1100005	Acid Phosphatase	Serum/Plasma		
1101000	Albumin	Serum/Plasma		
1103005	Alkaline Phosphatase	Serum/Plasma		
1105000	ALT/GPT	Serum/Plasma		
1107005	Amylase	Serum/Plasma/Urine		
1109000	AST/GOT	Serum/Plasma		
1110005	Bilirubin Direct	Serum		
1111010	Bilirubin Total	Serum		
1113000	Calcium Arsenazo	Serum/Plasma/Urine		
1115000	Calcium OCC	Serum/Plasma/Urine		
1116005	Chloride	Serum/Plasma/CSF		
1118005	Cholesterol	Serum/Plasma		
1119005	Cholinesterase	Serum/Plasma		
1120010	Creatinine Kinase BR	Serum/Plasma		
1121005	Creatinine Kinase-MB	Serum		
1123005	Creatinine	Serum/Plasma/Urine		
1126005	GGT	Serum/Plasma		
1129005	Glucose	Serum/Plasma/CSF/ Urine		
1133010	HDL-Cholesterol	Serum/Plasma		
1134015	Hemoglobin	WB		
1135105	Iron Cromazurol	Serum/Plasma		
1135005	Iron Ferrozine	Serum/Plasma		
1141010	LDH	Serum		
1133105	LDL-Cholesterol	Serum/Plasma		
1143010	Lipase	Serum/Plasma		
1144005	Magnesium	Serum/Plasma		
1148010	Phosphorus	Serum/Plasma/Urine		
1149005	Phosphorus UV	Serum/Plasma/Urine		
1150015	Potassium	Serum		
1162005	Protein (Urine/CSF)	Urine/CSF		
1151015	Sodium	Serum		
1137010	TIBC	Serum/Plasma		
1153005	Total Protein	Serum/Plasma		
1155005	Triglycerides	Serum/Plasma		

**Can be tested in KROMA, KROMA Plus, LIDA 500 and TEMIS.**



ALTAMEDICA INCORPORATED

Product List  
VETERINARY

*Dry Clinical Chemistry*

Code	Product Name	Description	Brand
N/A	Iris	Requires 140 uL of WB/Serum/Plasma; Release results after 9-12 mins.; With built-in printer and QR scanner	Linear

*Test Panel*

Code	Panel Name	Parameters	Packaging	Shelf-Life
N/A	General Panel	ALT, ALP, BUN, Crea, TBil, TP, Alb, Chol, Amy, Glu, Na, K, Na/K, Glob, A/G, URA/Crea		
N/A	Detailed Panel	ALT, AST, GGT, ALP, BUN, Crea, TBil, TP, Alb, Chol, Amy, Glu, AST/ALT, Glob, A/G, URA/Crea		
N/A	Pre-Operative Panel	ALT, ALP, BUN, Crea, TP, Alb, Glu, Na, K, Cl, Na/K, Glob, A/G, URA/Crea		
N/A	Renal Panel	BUN, Crea, Alb, Amy, Ca, Na, K, Cl, P, Na/K, Glob, A/G, URA/Crea		
N/A	Liber Panel	ALT, AST, GGT, ALP, BUN, TBil, TP, Alb, TBA, Chol, AST/ALT, Glob, A/G		
N/A	Electrolyte Panel	Ca, Na, K, Cl, P, Mg, Na/K		
N/A	Intensive Care Panel	ALT, BUN, Crea, Amy, LDH, Glu, CO <sub>2</sub> , Na, K, Cl, Na/K, URA/Crea		

*Coagulation*

Code	Product Name	Description	Parameters
3801100	CERES	Semi-automated coagulation machine with 2 test channels, 2 reagent positions and 4 sample positions	PT-HS, APTT, Fibrinogen
3801000	ARES	Semi-automated coagulation machine with 4 test channels, 6 reagent positions and 24 sample positions	PT-HS, APTT, Fibrinogen



ALTAMEDICA INCORPORATED

Product List  
VETERINARY

*Hematology*

Code	Product Name	Description	Brand
N/A	STEL3 VET	3-part differential hematology analyzer with 60 samples per hour throughput	Linear
N/A	STEL5 VET	5-part differential hematology analyzer with 60 samples per hour throughput	Linear

*Stel3 Vet and Stel5 Vet Reagents*

Code	Product Name	Packaging	Brand
N/A	Stel3 Vet Diluent		Linear Chemicals SLU
N/A	Stel3 Vet Lyse		Linear Chemicals SLU
N/A	Stel3 Vet Probe Cleanser		Linear Chemicals SLU
N/A	Stel5 Vet Diluent		Linear Chemicals SLU
N/A	Stel5 Vet LH Lyse		Linear Chemicals SLU
N/A	Stel5 Vet Diff Lyse		Linear Chemicals SLU
N/A	Stel5 Vet Probe Cleanser		Linear Chemicals SLU

*Electrolytes*

Code	Product Name	Description	Brand
N/A	AFT-500VET	Specialized electrolyte analyzer for veterinary; Test parameters Na, K, Cl, Ca, and pH; High speed test circle with less than 30 seconds per sample	Meizhou Cornley Hi-Tech Co., Ltd.
N/A	AFT-800VET	Specialized electrolyte analyzer for veterinary; Test parameters Na, K, Cl, Ca, and pH; High speed test circle with less than 25 seconds per sample	Meizhou Cornley Hi-Tech Co., Ltd.



ALTAMEDICA INCORPORATED

Product List  
VETERINARY

*AFT-500VET and AFT-800VET Reagents and Consumables*

Code	Product Name	Description	Packaging	Shelf-Life
10.10.000 501	AFT-500VET Calibration Standard Solution	Calibration Standard Solution for use with ISE (K <sup>+</sup> , Na <sup>+</sup> , Cl <sup>-</sup> , Ca <sup>++</sup> , pH) System		
10.10.162 004	Electrolytes Quality Control (Low/Middle/High) for AFT-500VET and AFT-800VET			
10.10.141 202	Conditioning Solution for AFT-500VET and AFT-800VET	Conditioning solution for use with Na <sup>+</sup> and pH electrodes for routine maintenance	5 pcs. per set	
10.10.151 302	Multi-enzyme De-Proteinizer for AFT-500VET and AFT-800VET	De-proteinizer for use with ISE measuring system for routine maintenance	3x2mL	
10.10.121 200	Refill Solution for ISE for AFT-500VET and AFT-800VET	Filling Solution for use only with K <sup>+</sup> , Na <sup>+</sup> , Cl <sup>-</sup> , Ca <sup>++</sup> , and pH Electrodes	5 pcs. per set	
10.10.131 201	Refill Solution Reference Electrode for AFT-500VET and AFT-800VET	Filling Solution for use only with Reference Electrode	20mL per bottle	
10.01.010 100	K Electrode for AFT-500VET		per piece	
10.01.020 100	Na Electrode for AFT-500VET		per piece	
10.01.030 100	Cl Electrode for AFT-500VET		per piece	
10.01.040 100	Ca Electrode for AFT-500VET		per piece	
10.01.060 100	pH Electrode for AFT-800VET		per piece	
10.01.070 100	Reference Electrode for AFT-500VET		per piece	



ALTAMEDICA INCORPORATED

Product List  
VETERINARY

Code	Product Name	Description	Packaging	Shelf-Life
10.10.000 014	AFT-800VET Calibration Standard Solution for ISE	Calibration Standard Solution for use with ISE (K+, Na+, Cl-, Ca++, pH) System	850mL	
10.10.950 092	AFT-800VET Calibration Standard Solution for ISE	Calibration Standard Solution for use with ISE (K+, Na+, Cl-, Ca++, pH) System	1250mL	
10.01.010 201	K Electrode for AFT-800VET		per piece	
10.01.020 200	Na Electrode for AFT-800VET		per piece	
10.01.030 101	Cl Electrode for AFT-800VET		per piece	1 year
10.01.040 201	Ca Electrode for AFT-800VET		3 x 2.0mL	1 year
10.01.060 200	pH Electrode for AFT-800VET		5 x 0.8mL	1 year
10.01.070 200	Reference Electrode for AFT-800VET			

*Disclaimer: All the products listed herein are subject to change from time to time without prior notice. . For inquiries , kindly contact customerservice@altamedica.com.ph.*